



Test kit, total chlorine, model CN-DT, 20-2000 mg/L

Productnr.: 2471100
EUR Prijs: Neem contact op
Verzenddatum n.t.b.

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Kit contains all required reagents and apparatus to perform analysis in the field

Specificaties

Aantal testen:	100
Koffer design:	F
Meetbereik:	20 - 2,000 mg/L Cl ₂
Methode:	Digital Titrator/Thiosulfate
Opmerking:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Parameter:	Chlorine, mid and high range - As total Cl ₂
Platform:	Digitale Titrator
Smallest increment:	0.2 - 5.0
Testkit model:	CN-DT
Verzend gewicht:	4